

ÉCOLE POLYTECHNIQUE FÉDÉRALE DE LAUSANNE

MASTER THESIS

Evaluation of optical aberrations using Phase Diversity

Author:
Jordan VOIRIN

Supervisors:
Dr. Laurent JOLISSAINT
Dr. Jean-Paul KNEIB

*A thesis submitted in fulfillment of the requirements
for the degree of Master in Applied Physics*

in the

Astrophysics laboratory
Basic Sciences

November 1, 2017

Declaration of Authorship

I, Jordan VOIRIN, declare that this thesis titled, "Evaluation of optical aberrations using Phase Diversity" and the work presented in it are my own. I confirm that:

- This work was done wholly or mainly while in candidature for a research degree at this University.
- Where any part of this thesis has previously been submitted for a degree or any other qualification at this University or any other institution, this has been clearly stated.
- Where I have consulted the published work of others, this is always clearly attributed.
- Where I have quoted from the work of others, the source is always given. With the exception of such quotations, this thesis is entirely my own work.
- I have acknowledged all main sources of help.
- Where the thesis is based on work done by myself jointly with others, I have made clear exactly what was done by others and what I have contributed myself.

Signed:

Date:

“Thanks to my solid academic training, today I can write hundreds of words on virtually any topic without possessing a shred of information, which is how I got a good job in journalism.”

Dave Barry

ÉCOLE POLYTECHNIQUE FÉDÉRALE DE LAUSANNE

Abstract

Physics
Basic Sciences

Master in Applied Physics

Evaluation of optical aberrations using Phase Diversity

by Jordan VOIRIN

The Thesis Abstract is written here (and usually kept to just this page). The page is kept centered vertically so can expand into the blank space above the title too...

Acknowledgements

The acknowledgments and the people to thank go here, don't forget to include your project advisor...

Contents

Declaration of Authorship	iii
Abstract	vii
Acknowledgements	ix
1 Introduction	1
2 Phase Diversity Experiment	3
2.1 Theoretical Background	3
2.2 Experimental Setup	3
2.2.1 Light source	3
2.2.2 Entrance pupil	4
2.2.3 Pupil imaging system	4
2.2.4 Acquisitions	4
2.3 Results	5
2.3.1 Astigmatism	5
A Optical Component Datasheets	7
A.1 Pigtailed laser diode	8
A.2 Converging lens A220TM-A, $f = 11$ mm	9
A.3 Pinhole $10\ \mu\text{m}$	10
A.4 Converging lens AL100200, $f = 200$ mm	11
A.5 Converging lens AC254-100-A, $f = 100$ mm	12
A.6 Ximea Camera, MQ013MG-E2	13
A.7 Shack-Hartmann wavefront sensor, WFS150-5C	14
Bibliography	15

List of Figures

2.1	Schema of the experimental setup	3
2.3	Source schema and pinhole effect on the beam.	4
2.2	Wavefront curvature	4

List of Tables

2.1 Optical Components	3
----------------------------------	---

List of Abbreviations

WFS WaveFront Sensor

Physical Constants

Speed of Light $c_0 = 2.997\,924\,58 \times 10^8 \text{ m s}^{-1}$ (exact)

List of Symbols

a	distance	m
P	power	W (J s ⁻¹)
ω	angular frequency	rad

For/Dedicated to/To my...

Chapter 1

Introduction

???

Chapter 2

Phase Diversity Experiment

2.1 Theoretical Background

2.2 Experimental Setup

The design of the experiment was already done by Bouxin (2017). The system is built according to her plans and specifications.

The experiment is mounted on a pressurized legs optical table. The assembly contains six main components : a light source, an entrance pupil, an imaging system, a converging lens to focus the beam on the camera, a camera and a wavefront sensor.

2.2.1 Light source

The final application of the phase diversity will be to characterize the optical aberrations induced by the imperfect optical path to a scientific detector of a telescope. For this reason, the light source has to simulate a distant




FIGURE 2.1: Experimental setup schema with the relevant distances, (Bouxin, 2017).

TABLE 2.1: Optical Components

#	Components	Model	Reference
1	Pigtailed laser diode	Thorlabs, LPS-635-FC	A.1
2	Converging lens, $f = 11$ mm	Thorlabs, A220TM-A	A.2
3	Pinhole, $10 \mu\text{m}$	Thorlabs, P10S	A.3
4	Converging lens, $f = 200$ mm	Thorlabs, AL100200	A.4
5	3.2 mm Hole milled in metal sheet
6	Converging lens, $f = 100$ mm	Thorlabs, AC254-100-A	A.5
7	Converging lens, $f = 80$ mm		
8	Camera CMOS	Ximea, MQ013MG-E2	A.6
9	Converging lens, $f = 100$ mm		
10	Shack-Hartman WFS	Thorlabs, WFS150-5C	A.7





FIGURE 2.3: Source schema and pinhole effect on the beam.

star aberration-free wavefront. A distant star wavefront is considered planar since the object distance, z , is far greater than the telescope size, r , see Fig. 2.2. The source of our experiment must then be characterized by planar wavefront.

In order to obtain such a planar wavefront at the entrance pupil, the light source consist of a "pigtailed laser diode", a $f=11\text{mm}$ converging lens, a pinhole and a $f=200\text{ mm}$ converging lens, see Table 2.1. The pigtailed laser diode emits a Gaussian beam centred at 637.5 nm slightly diverging. The converging lens concentrates the beam at the center of the $10\mu\text{m}$ pinhole to filter the noise. The second converging lens collimates the beam, obtaining a collimated beam with a planar wavefront, see Fig. 2.3a and 2.3b.

FIGURE 2.2: Wavefront curvature for different source's distances, z . r represents the characteristic size of the arc of interest.

2.2.2 Entrance pupil

The entrance pupil of our optical system is a circular aperture of 3.2 mm diameter placed after the collimating lens of the light source. It is milled in a metal plate and centred in his support, to avoid positioning with a XY table. The diameter is chosen in available material to fit in the different detector's surfaces.

2.2.3 Pupil imaging system

The phase diversity technique requires PSFs images as input, which means that the beam as to be focused onto the detector surface. To analyse the aberration in the pupil plane, one needs to focus an image of the beam passing through the entrance pupil. The simplest assembly to achieve this goal is the 4F system, which consist of two converging lenses of focal 100 mm . The two lenses are separated by 200 mm , see Fig. 2.1. This places the image of the entrance pupil 100 mm after the second converging lens.

2.2.4 Acquisitions

The image of the entrance pupil, obtained with the 4F system, is focused onto a CMOS Ximea camera by a $f = 80\text{ mm}$ converging lens to acquire the PSFs for the phase diversity wavefront retrieval. The camera has surface composed by 1280×1024

pixels of $5.3\ \mu\text{m}$, see Appendix A.6. It is mounted on sliding support in order to be able to acquire in/out-of-focus images. A beam splitter is placed in the converging beam to separate it in two. The second beam is collimated and a Shack-Hartman WFS is placed on the entrance pupil image plane, to check the results of the phase diversity wavefront retrieval. The Shack-Hartman WFS has a 39×31 lenslet grid and a CCD with a resolution of 1280×1024 pixels of $4.65\ \mu\text{m}$, see Appendix A.7.

2.3 Results

This section presents the results of the phase diversity experiment, with the introduction of different sources of aberration.

2.3.1 Parallel plane plate

The first source of aberration studied in this work is a tilted parallel plane plate which is used as a calibrated source of astigmatism.


Appendix A

Optical Component Datasheets

A.1 Pigtailed laser diode

THORLABS

Pigtailed Laser Diode, SMF



Description

Thorlabs' Single Mode Pigtailed Laser Diodes are standard TO-packaged diodes that have been pigtailed to a 1 m long single mode fiber with an FC/PC connector. Each unit is tested before shipment. Please refer to the unit-specific test datasheet for optimal operating parameters.

Specifications


LPS-635-FC Specifications	
LD Reverse Voltage (Max)	2 V
PD Reverse Voltage (Max)	30 V
Optical Output Power	2.5 mW (Typ.) 3.5 mW (Max)
Operating Temperature	0 to 50 °C
Storage Temperature	-10 to 65 °C
Pin Code	9A
Laser Diode	HL6320G
Fiber	SM600
Connector	FC/PC

	Min	Typ.	Max
Wavelength	625 nm	635 nm	640 nm
Threshold Current*	20 mA	50 mA	75 mA
Slope Efficiency*	0.13 mW/mA	0.18 mW/mA	-
Operating Current @ $P_O = 2.5 \text{ mW}^*$	-	70 mA	95 mA
Operating Voltage @ $P_O = 2.5 \text{ mW}^*$	-	2.2 V	2.7 V
Monitor Current @ $P_O = 2.5 \text{ mW}^*$	0.05 mA	0.17 mA	0.3 mA


*Temperature = 25 °C

Drawing


Pigtail Bottom View




Side View




Diode Bottom View



Pin Diagram





NOTICE

To avoid equipment damage from electrostatic discharge: Wear ESD wriststrap when handling this device.


US, Canada, & South America: +1-973-300-3000 | France: +33 (0) 970 444 844 | Europe: +49 (0) 8131-5956-0 | UK & Ireland: +44 (0)1353-654440
 Brazil: +55-16-3413 7062 | Scandinavia: +46-31-733-30-00 | Japan & Asia: +81-3-5979-8889 | China: +86 (0)21-60561122

June 6, 2013
5175-S01, Rev B

www.thorlabs.com


Source : www.thorlabs.com

A.2 Converging lens A220TM-A, f = 11 mm




Source : www.thorlabs.com

A.3 Pinhole 10 μm




Source : www.thorlabs.com

A.4 Converging lens AL100200, f = 200 mm



Source : www.thorlabs.com

A.5 Converging lens AC254-100-A, $f = 100$ mm



Source : www.thorlabs.com

A.6 Ximea Camera, MQ013MG-E2



Specifications:

Resolution:	1.3 MP 1280 × 1024 pixels
Sensor type:	CMOS Matrix B/W
Sensor model:	e2V EV76C560 ABT-EQV
Sensor size:	1/1.8"
Sensor active area:	6.9 × 5.5 mm
Pixel size:	5.3 μm
Bits per pixel:	8, 10
Dynamic range:	60 dB
Frame rates:	60 fps
On-chip binning:	1x1, 2x2
Image data interface:	USB 3.0
Data I/O:	GPIO IN, OUT
Power requirements:	0.9 Watt
Lens mount:	C or CS Mount
Weight:	26 grams
Dimensions WxHxD:	26 x 26 x 26 mm
Operating environment:	50 °C
Customs tariff code:	8525.80 30 (EU) / 8525.80 40 (USA)
ECCN:	EAR99

Source : www.ximea.com/en/products/usb3-vision-cameras-xiq-line/mq013mg-e2

A.7 Shack-Hartmann wavefront sensor, WFS150-5C

8 Appendix

8.1 Technical Data

8.1.1 WFS150/300

Item #	WFS150-5C	WFS150-7AR	WFS300-14AR
Microlenses			
Microlens Array	MLA150M-5C	MLA150M-7AR	MLA300M-14AR
Substrate Material	Fused Silica (Quartz)		
Number of Active Lenslets	Software Selectable		
Max. Number of Lenslets	39 x 31		19 x 15
Camera			
Sensor Type	CCD		
Resolution	max. 1280 x 1024 pixels, Software Selectable		
Aperture Size	5.95 mm x 4.76 mm		
Pixel Size	4.65 μm x 4.65 μm		
Shutter	Global		
Exposure Range	79 μs - 65 ms		
Frame Rate	max. 15 Hz		
Image Digitization	8 bit		
Wavefront Measurement			
Wavefront Accuracy ¹⁾	$\lambda/15$ rms @ 633 nm		$\lambda/50$ rms @ 633 nm
Wavefront Sensitivity ²⁾	$\lambda/50$ rms @ 633 nm		$\lambda/150$ rms @ 633 nm
Wavefront Dynamic Range ³⁾	> 100 λ @ 633 nm		> 50 λ @ 633 nm
Local Wavefront Curvature ⁴⁾	> 7.4 mm	> 10.0 mm	> 40.0 mm
External Trigger Input			
Save Static Voltage level	0 to 30 V DC		
LOW Level	0.0 V to 2.0 V		
HIGH Level	5.0 V to 24 V		
Input current	> 10 mA		
Min Pulse Width	100 μs		
Min. Slew Rate	35 V / msec		
Common Specifications			
Optical Input	C-Mount		
Power Supply	<1.5 W, via USB		
Operating Temperature Range ⁵⁾	+5 to +35 $^{\circ}\text{C}$		
Storage Temperature Range	-40 to 70 $^{\circ}\text{C}$		
Warm-Up Time for Rated Accuracy	15 min		
Dimensions (W x H x D)	32.0 mm x 40.4 mm x 45.5 mm		
Weight	0.1 kg		

¹⁾ Absolute accuracy using internal reference. Measured for spherical wavefronts of known RoC.

²⁾ Typical relative accuracy. Achievable after, and with respect to a user calibration, 10 image averages

³⁾ Over entire aperture of wavefront sensor

⁴⁾ Radius of wavefront curvature over single lenslet aperture

⁵⁾ non-condensing

All technical data are valid at 23 \pm 5°C and 45 \pm 15% rel. humidity (non condensing)

Bibliography

Bouxin, A. (2017). "Phasor diversity to measure the static aberrations of an optical system". MA thesis. HEIG-VD.

Thorlabs (2017). *Principles of Spatial Filters*. Thorlabs. URL: https://www.thorlabs.com/newgrouppage9.cfm?objectgroup_id=1400.